



SWEAT STOPS HERE™

PERMANENTLY REDUCE UNDERARM
SWEAT, ODOR AND HAIR^{1,2}

CONSULTATION GUIDE




EVERYONE SWEATS

1

For some, sweating is excessive


For others, sweating and odor are simply a bother – managed daily with antiperspirant and deodorant



SWEAT FACTS

ANATOMY OF SWEAT AND ODOR GLANDS

2



DO I NEED UNDERARM SWEAT?

While your body does need sweat glands to cool itself, **your body contains 2-4 million sweat glands.** Only about 2% are found in your underarms.


Eliminating this 2% of sweat glands does not affect your body's ability to cool itself.



HOW MIRADRY WORKS

4

miraDry delivers targeted thermal energy to safely eliminate sweat, odor, and hair.^{1,2}



A COMFORTABLE 1-HOUR TREATMENT

5

1. Numb treatment area



2. Apply temporary tattoo



3. Treat with miraDry system



82%

Average reduction in underarm sweat³

MIRADRY IS FOR EVERYONE

Social Situations	Avoid sweating during social events
Professionals	Avoid sweating during meetings
Natural Lifestyle	Avoid chemicals/parabens in antiperspirants
Wardrobe Conscious	Tired of yellowed shirts or buying new clothes
Aesthetic Minded	Enjoys taking care of him/herself
Young Adults	Sweat affects his/her confidence
Solution Seeker	Sweats more than average



WHO IS MIRADRY FOR?

WHAT PATIENTS CAN EXPECT

7

Comfortable, one-hour procedure

Common side effects

- 💧 Swelling
- 💧 Soreness
- 💧 Altered sensation in or around underarm

Little to no downtime

- 💧 Use ice and ibuprofen to reduce swelling
- 💧 Wait a few days before resuming rigorous exercise and activity



WHAT TO EXPECT

THE BENEFITS OF MIRADRY ARE CLEAR

8

Permanently reduces underarm sweat, odor and hair

Immediate results in as little as one treatment

Non-surgical

Little to no downtime

FDA-cleared




patient
satisfaction³



PROVEN RESULTS

Before miraDry



After miraDry



Photos courtesy of SkinSpeaks, Edina, MN

Actual miraDry patient who took photos after a workout before and after her miraDry treatment.

As a TV host and former Miss America Runner Up, wardrobe changes are a big part of my life. With miraDry, I don't have to ever worry about unsightly sweat stains on camera. What a relief!

Crystal Lee

TV Host

Photograph by Al Ponce



IMPORTANT SAFETY INFORMATION

The miraDry treatment is not for you if you have or need the following:

- Heart pacemakers and other electronic device implants.
- Supplemental oxygen.
- Known resistance to or history of intolerance of local anesthesia including lidocaine and epinephrine.

Use caution before undergoing the miraDry treatment under the following conditions, the effects of which have not been studied:

- Pregnancy
- Underlying skin conditions
- Previous axillary surgery

COMMON AND RARE SIDE EFFECTS

The following effects can occur during and after a miraDry treatment. In most cases, these side effects will gradually go away. In rare cases, it can last for several months.

DURING A TREATMENT

- Due to the anesthesia, you may experience bruising at the injection site. You may also experience shaking, numbness or tingling in the arm, lasting less than 24 hours.

IMMEDIATELY AFTER A TREATMENT

- Swelling, redness, temporary altered sensation, tingling, soreness, weakness, tight banding, pain, or bumps under the skin in the treated area and/or upper arm.
- Discomfort, tenderness or pain in the underarm is typically treatable with non-prescription medications such as ibuprofen. In rare cases, prescription medications may be needed.

RARE SIDE EFFECTS

- Hyperpigmentation (darkening of skin), burns, skin infections, rashes, and altered sweating in other areas of the body may occur.

RESULTS:

- Results vary from person to person. You may decide that additional treatments are necessary to achieve your desired outcome. Although highly unlikely, it is possible that you will not experience any noticeable result from the procedure.

Learn more at miraDry.com

References

1. Li et al. Antigen expression of human eccrine sweat glands. J Cutan Pathol 36: 318-324 (2009). Lupin et al. Long-term efficacy and quality of life assessment for treatment of axillary hyperhidrosis with a microwave device. Dermatol Surg 40: 805-807 (2014).
2. Data on file, miraDry, Inc.
3. Hong et al. Clinical Evaluation of a Microwave Device for Treating Axillary Hyperhidrosis. Dermatol Surg 38:728-735 (2012).

miraDry, Inc. - 2790 Walsh Avenue – Santa Clara, CA 95051

Patient results and experience may vary. The miraDry system is FDA-cleared for the treatment of unwanted underarm sweat and odor, as well as the permanent removal of unwanted underarm hair. miraDry® and Sweat Stops Here™ are trademarks of miraDry, Inc. © 2018 miraDry, Inc. All rights reserved.

MK0696.A 07/18

CONSULTATIONS AVAILABLE

miraDry®
The sweat stops here.™